

Creating modified atmospheres



Anaerobic

Don Whitley Scientific has designed and developed anaerobic workstations continuously since 1980.

We have one of the most comprehensive product ranges for routine and research applications.

Microaerobic

In 1983, Don Whitley Scientific developed the world's first variable atmosphere workstation for microbiological applications.

Hypoxic

Designed specifically for tissue and cell culture applications, our hypoxic workstations provide precise and accurate control of oxygen, carbon dioxide, temperature and humidity levels.

Jars

Our high quality jars incorporate a patented valve system and can be fitted with a low temperature catalyst sachet. Our jars can be used with gas generating envelopes or prepared using the Whitley Jar Gassing System.

Jar Gassing Systems


The Whitley Jar Gassing System creates microaerobic or anaerobic conditions in seconds for just a few pence, without the need for gas generation kits and without creating chemical waste.

Workstations

Samples can be introduced, manipulated, incubated, examined and removed without loss of controlled environmental conditions.

Isolation rates are significantly increased when samples are placed in a workstation immediately after collection and remain there until identification.*

Anaerobic workstations/1



Features	Whitley 500	Whitley 500 + Airlock	Whitley 1000	Whitley 1000 + Airlock	Features
Chamber Capacity	540 / 700	540 / 700	1080 / 1400	1080 / 1400	Chamber Capacity
Port / Airlock Capacity	20 plates per port (40)	20 plates per port + 90 plates via airlock (130)	20 plates via ports (80)	20 plates per port + 90 plates via airlock (170)	Port / Airlock Capacity
Porthole System	Motorised	Motorised	Motorised	Motorised	Porthole System
Gas Supplies	ANO ₂ / N ₂	H ₂ / CO ₂ / N ₂	ANO ₂ / N ₂	H ₂ / CO ₂ / N ₂ or ANO ₂ / N ₂	Gas Supplies
Footswitch	Wireless	Wireless	Wireless	Wireless	Footswitch
Auto Sleeve Gassing	-	-	-	-	Auto Sleeve Gassing
Internal Mains Socket	○	○	○	○	Internal Mains Socket
Storage Trays	-	-	-	-	Storage Trays
Lighting	●	●	●	●	Lighting
Inspection Lamp	○	○	○	○	Inspection Lamp
Single Sample Entry	○	○	○	○	Single Sample Entry
ANO ₂ Indicator	○	○	○	○	ANO ₂ Indicator
ANO ₂ Conditions Monitor	-	-	-	-	ANO ₂ Conditions Monitor
Refrigeration	○	○	○	○	Refrigeration
Data Logging	-	-	-	-	Data Logging
Airlock Cycle Time	-	6 minutes	-	6 minutes	Airlock Cycle Time
Extra Cable Glands	○	○	○	○	Extra Cable Glands
Microscope	○	○	○	○	Microscope
Automatic Dehumidifier	●	●	●	●	Automatic Dehumidifier
Automatic Humidifier	-	-	-	-	Automatic Humidifier
Removable Front	-	-	-	-	Removable Front
Workstation Trolley	○	○	○	○	Workstation Trolley
Dimensions (w/d/h)	1040 / 760 / 800	1570 / 760 / 840	1885 / 760 / 800	2415 / 760 / 840	Dimensions (w/d/h)
Weight (lbs/kg)	220 / 100	340 / 155	385 / 175	506 / 230	Weight (lbs/kg)

Notes : ● Fitted as standard ○ Option available - Not applicable

*Capacity based upon using 90mm Petri dishes. Shown as working / short term maximum with limited working space

Anaerobic workstations/2



Features	DG250	A35	A45	A85	Features
Chamber Capacity	270 / 400	400 / 600	600 / 900	540 / 700	Chamber Capacity
Port / Airlock Capacity	10 plates per port (20)	40 plates via airlock	40 plates via airlock	20 plates per port + 90 plates via airlock (130)	Port / Airlock Capacity
Porthole System	Manual	Instant / Manual	Instant Access Porthole System	Motorised	Porthole System
Gas Supplies	ANO ₂ / N ₂	ANO ₂ / N ₂	ANO ₂ / N ₂	ANO ₂ / N ₂	Gas Supplies
Footswitch	Hardwired	Wireless	-	Wireless	Footswitch
Auto Sleeve Gassing	○	○	○	○	Auto Sleeve Gassing
Internal Mains Socket	○	○	○	○	Internal Mains Socket
Storage Trays	○	○	○	-	Storage Trays
Lighting	○	●	●	●	Lighting
Inspection Lamp	●	○	○	○	Inspection Lamp
Single Sample Entry	○	○	○	○	Single Sample Entry
ANO ₂ Indicator	○	-	-	○	ANO ₂ Indicator
ANO ₂ Conditions Monitor	-	○	○	○	ANO ₂ Conditions Monitor
Refrigeration	-	-	-	○	Refrigeration
Data Logging	-	○	○	○	Data Logging
Airlock Cycle Time	-	60 seconds	60 seconds	TBA	Airlock Cycle Time
Extra Cable Glands	○	○	○	○	Extra Cable Glands
Microscope	○	○	○	○	Microscope
Automatic Dehumidifier	●	●	●	●	Automatic Dehumidifier
Automatic Humidifier	-	○	○	-	Automatic Humidifier
Removable Front	Full lift off lid	○	○	-	Removable Front
Workstation Trolley	○	○	○	○	Workstation Trolley
Dimensions (w/d/h)	810 / 760 / 635	1255 / 720 / 710	1650 / 720 / 710	1570 / 760 / 840	Dimensions (w/d/h)
Weight (lbs/kg)	141 / 64	220 / 100	TBC	330 / 150	Weight (lbs/kg)

Notes : ● Fitted as standard ○ Option available - Not applicable *Capacity based upon using 90mm Petri dishes. Shown as working / short term maximum with limited working space

Microaerobic workstations



Features	Whitley VA500	Whitley VA500 + Airlock	Whitley VA1000	Whitley VA1000 + Airlock	Features
Chamber Capacity	540 / 700	540 / 700	1080 / 1400	1080 / 1400	Chamber Capacity
Port / Airlock Capacity	20 plates per port (40)	20 plates per port + 90 plates via airlock (130)	20 plates per port (80)	20 plates per port + 90 plates via airlock (170)	Port / Airlock Capacity
Porthole System	Motorised	Motorised	Motorised	Motorised	Porthole System
Gas Supplies	Up to 4 separate cylinders of H ₂ /N ₂ , CO ₂ , N ₂ , Air				Gas Supplies
Footswitch	Wireless	Wireless	Wireless	Wireless	Footswitch
Auto Sleeve Gassing	-	-	-	-	Auto Sleeve Gassing
Internal Mains Socket	○	○	○	○	Internal Mains Socket
Storage Trays	-	-	-	-	Storage Trays
Lighting	●	●	●	●	Lighting
Inspection Lamp	○	○	○	○	Inspection Lamp
Single Sample Entry	○	○	○	○	Single Sample Entry
ANO ₂ Indicator	-	-	-	-	ANO ₂ Indicator
ANO ₂ Conditions Monitor	-	-	-	-	ANO ₂ Conditions Monitor
Refrigeration	○	○	○	○	Refrigeration
Data Logging	-	-	-	-	Data Logging
Airlock Cycle Time	-	5 minutes	-	5 minutes	Airlock Cycle Time
Extra Cable Glands	○	○	○	○	Extra Cable Glands
Microscope	○	○	○	○	Microscope
Automatic Dehumidifier	●	●	●	●	Automatic Dehumidifier
Automatic Humidifier	-	-	-	-	Automatic Humidifier
Removable Front	-	-	-	-	Removable Front
Workstation Trolley	○	○	○	○	Workstation Trolley
Dimensions (w/d/h)	1040 / 760 / 800	1570 / 760 / 840	1885 / 760 / 800	2415 / 760 / 840	Dimensions (w/d/h)
Weight (lbs/kg)	242 / 110	363 / 165	407 / 185	528 / 240	Weight (lbs/kg)

● Fitted as standard ○ Option available - Not applicable

*Capacity based upon using 90mm Petri dishes. Shown as working / short term maximum with limited working space

Hypoxic workstations



Features	H35 Hypoxystation	H85 Hypoxystation	Features
O ₂ Profiling	<input type="radio"/>	<input type="radio"/>	O ₂ Profiling
CO ₂ Monitoring	<input type="radio"/>	<input type="radio"/>	CO ₂ Monitoring
Data Logging	<input type="radio"/>	<input type="radio"/>	Data Logging
Automatic Dehumidifier	<input checked="" type="radio"/>	<input checked="" type="radio"/>	Automatic Dehumidifier
Automatic Humidifier	<input type="radio"/>	<input type="radio"/>	Automatic Humidifier
Removable Front	<input type="radio"/>	-	Removable Front
Chamber Capacity	326 litres	300 litres	Chamber Capacity
Airlock Capacity	12 litres	36 litres	Airlock Capacity
Gas Supplies	CO ₂ / Air / N ₂	CO ₂ / Air / N ₂	Gas Supplies
Porthole System	Conventional Manual or Instant Access Porthole System	Motorised	Porthole System
Footswitch	Wireless	Wireless	Footswitch
Auto Sleeve Gassing	<input type="radio"/>	<input type="radio"/>	Auto Sleeve Gassing
Internal Mains Socket	<input type="radio"/>	<input type="radio"/>	Internal Mains Socket
Storage Trays	<input type="radio"/>	-	Storage Trays
Lighting	<input checked="" type="radio"/>	<input checked="" type="radio"/>	Lighting
Inspection Lamp	<input type="radio"/>	<input type="radio"/>	Inspection Lamp
Single Sample Entry	<input type="radio"/>	<input type="radio"/>	Single Sample Entry
Vacuum Take-off	<input type="radio"/>	<input type="radio"/>	Vacuum Take-off
Airlock Cycle Time	60 seconds	3.5 minutes	Airlock Cycle Time
Extra Cable Glands	<input type="radio"/>	<input type="radio"/>	Extra Cable Glands
Microscope	<input type="radio"/>	<input type="radio"/>	Microscope
Workstation Trolley	<input type="radio"/>	<input type="radio"/>	Workstation Trolley
Dimensions (w/d/h)	1255 / 720 / 710	1570 / 760 / 840	Dimensions (w/d/h)
Weight (lbs/kg)	231 / 105	330 / 150	Weight (lbs/kg)

● Fitted as standard

○ Option available

- Not applicable

Maintaining good anaerobic conditions

Stable Temperature

Temperature distribution within the workstation should be uniform. DWS workstations maintain temperature within $\pm 0.6^{\circ}\text{C}$ of the set point across the working and incubation areas of the chamber.

No Oxygen

Anaerobic gas mixture contains 10% hydrogen, 10% carbon dioxide and 80% nitrogen. Residual oxygen is removed by combining with hydrogen in the presence of a catalyst to produce water.



Catalyst + Anotox™

Catalyst and Anotox must be in good condition and not allowed to get wet. The typical life expectancy of catalyst is between 6-12 months depending on the types of bacteria being grown. Anotox is essential for maintaining ideal internal conditions by absorbing volatile fatty acids and hydrogen sulphide, both of which can be byproducts of the growth of anaerobic bacteria. Anotox has been shown to substantially prolong catalyst life. Don Whitley Scientific recommends that catalyst and Anotox are replaced together.

Relative Humidity Control

To maintain ideal bacterial growth conditions, agar must not be allowed to dry out during sample incubation. Excessive internal moisture may promote bacterial cross contamination.

A Sealed Workstation Chamber

Workstations are operated at a slight positive pressure so that any small leaks allow anaerobic gas to leak out, instead of allowing oxygen to enter the workstation.

Maintaining good hypoxic conditions

A Sealed Workstation Chamber


Workstations are operated at a slight positive pressure so that any small leaks allow internal gas to leak out, instead of allowing air to enter the workstation.

Stable Temperature

Temperature distribution within the workstation should be uniform. DWS workstations maintain temperature within $\pm 0.6^{\circ}\text{C}$ of the set point across the working and incubation areas of the chamber.

Oxygen Level Control

An automatic feedback system controls oxygen in 0.1% increments from 0.1% to 20%. The oxygen cell is calibrated in situ. The maximum oxygen concentration possible is directly related to the percentage of carbon dioxide selected. There must be a minimum of 80% nitrogen in the system with the balance available to split between oxygen and carbon dioxide.



Relative Humidity Control

An automatic dehumidification system is fitted as standard. Moisture may need to be added to the atmosphere to maintain required levels of relative humidity. An optional humidifier can be specified.

Carbon Dioxide Level Control

An optional automatic feedback system controls carbon dioxide in 0.1% increments from 0.1% to 15%. The maximum carbon dioxide concentration possible is directly related to the percentage of oxygen selected. There must be a minimum of 80% nitrogen in the system with the balance available to split between carbon dioxide and oxygen.